



USER MANUAL



WE HAVE YOU COVERED

Our passionate team brings you innovative, quality mobility solutions. Our products are reliable, sleek, and designed with simplicity so that you can rest assured that they work when you need them.

Thanks to our out of the box thinking, we are vigilant to meet your mobility needs with our simple and practical solutions. We recognize the important role that safety plays in consumer choices and are proud to offer a line of products that maintain the safety features found in your vehicle. Backed by third party verification for compliance to applicable Federal Motor Vehicle Safety Standards and Canadian Motor Vehicle Safety Standards, our products are engineered for safety and durability.

Our commitment to quality and safety along with our genuine desire to meet your mobility needs makes for a product we are sure you will be proud to own.





ASENTTO

USER GUIDE

To get into your vehicle

Your new **ASENTO** assists your transfer into your vehicle.

Follow the next few steps:

- Open the door and fold down the transfer board.
- Position yourself and your wheelchair so as to transfer as safely as possible onto the **ASENTO**. Be sure to put the brakes on your wheelchair to prevent it from moving.
- Once you have transferred onto the **ASENTO**, make sure you lower the security arm in place, around you. This security arm was designed to help you keep your balance while the **ASENTO** is ascending or descending.



USER GUIDE (CONTINUED)

To get into your vehicle

- Pick up the hand held pendant, press and hold the **UP** button so that the **ASENTO** brings you up to the height of your vehicle's seat.
- Once the **ASENTO** has reached the desired height that will offer you the best position to transfer, release the **UP** button. You can then transfer over into the vehicle. At any time during the motions you can release the **UP** button and the **ASENTO** will stop, allowing you to reposition. Once you are ready to continue, just press and hold the **UP** button to resume your ascent.
- After completing your transfer, press and hold the **DOWN** button to lower the **ASENTO** back to the floor of your vehicle. Once the **ASENTO** has reached its lowest point, release the button and stow the security arm first, followed by the transfer board, by folding it upwards. You can now close your vehicle's door.



USER GUIDE (CONTINUED)

To get out of your vehicle

Follow the next few steps:

- Remove your seat belt and open the door.
- Fold down the transfer board and the security arm.
- Pick up the hand held pendant.
- Press and hold the **UP** button until the **ASENTO** has reached the desired height for you to transfer onto it and then release the button.
- Transfer onto the **ASENTO** and use the security arm for balance. Press and hold the **DOWN** button until the **ASENTO** has reached the desired height that will offer you the best position to transfer and release the **DOWN** button.
- Stow the hand held pendant.
- Stow the security arm before attempting to transfer off the **ASENTO**.
- Transfer off the **ASENTO** and fold the transfer board back up in place and then close the vehicle's door.



WARNINGS

- WHEN CLOSING THE DOOR, BE SURE THAT THE TRANSFER BOARD IS FOLDED UP IN POSITION.
- THE TRANSFER BOARD MUST BE AT ITS LOWEST POINT POSSIBLE AND FOLDED UP BEFORE CLOSING THE DOOR.
- NEVER TRY TO OPERATE THE **ASENTO** WHILE THE VEHICLE DOOR IS CLOSED OR IF THE VEHICLE IS IN MOTION.
- **⚠WARNING:** This product can expose you to chemicals including benzene, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

POWER LOSS

If the 12 volt power supply is interrupted or removed, the **ASENTO** will return to programming mode. When the **ASENTO** is in programming mode it will be in reduced power, and will therefore beep and move very slowly. To reprogram the **ASENTO**, you must not be sitting on the unit. Simply run the transfer board all the way down until it stops and quickly beeps three (3) times. Then, run the board all the way up until it stops and quickly beeps three (3) times. The **ASENTO** will return to full power. Nothing further is required.

For more details, watch the video here:



3 YEAR LIMITED WARRANTY

ASENTO

Adapt Solutions # 1 Ltd. (Adapt Solutions), warrants to the original purchaser of an **ASENTO** that the equipment is free from defects in material and workmanship for a period of three years from date of purchase.

During the first year of the warranty, Adapt Solutions will supply the replacement parts as well as a pre-set monetary amount (determined by Adapt Solutions) for the repair if a defect in materiel or workmanship is discovered. After the initial year of this warranty, only parts and components are covered. This warranty does not cover labour and other services after the initial year. Freight and other related repair charges will be the responsibility of the original purchaser.

The only remedy for a defect in one of Adapt Solutions products (**LINK; XL-BASE; XL-SEAT; HI-LIFT; SPEEDY-LIFT; XL-BOARD; POWER-PULL; ASENTO**) shall be the repair or the replacement, at the discretion of Adapt Solutions, of the defective part or component. If repair or replacement is not commercially practical or cannot be timely made, Adapt Solutions may decide to refund the purchase price of the equipment instead of repairing or replacing the original equipment.

In no event shall Adapt Solutions be responsible for indirect, incidental or consequential damages, whether such damages arise from claims based on contract, warranty, tort (including negligence), strict liability or product liability.

All implied warranties, including any warranty of merchantability or fitness for a particular purpose, are limited in their duration to the length of the warranty stated above for the affected component.

This warranty is to the original purchaser only, and excludes product damage due to installation error, product misuse, product abuse, accidents, physical damage, damage in shipment, modifications not made by Adapt Solutions, or repairs undertaken by anyone other than authorized distributors.

TO OBTAIN WARRANTY PARTS OR REIMBURSEMENT

You must follow these procedures:

1. Obtain warranty authorization by calling your local Adapt Solutions dealer or Adapt Solutions directly at 1-866-641-0419.
2. Return the faulty Adapt Solutions component/equipment, freight prepaid, to the address provided by your Adapt Solutions dealer or Adapt Solutions with proof of purchase indicating the date purchased.

Adapt Solutions will pay for shipping back to the purchaser within the continental United States and Canada if a defect in materiel or workmanship is discovered. Return freight and repair charges will be the responsibility of the purchaser if the problem is not covered by warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state/ province to province. Adapt Solutions specifically does not authorize any person to extend the time or scope of this warranty.

For further information regarding this limited warranty, please contact us by calling 1-866-641-0419 or by writing at the following address:

**Adapt Solutions / Warranty Department
145 Damase-Breton
St-Lambert-de-Lauzon, Québec
G0S 2W0 Canada**

WARRANTY REGISTRATION

Adapt Solutions #1 Ltd located in St-Lambert-de-Lauzon, Québec, Canada, warrants this product to be free from any manufacturer's defect and workmanship for a period of three years from the date of purchase. This product has to be installed by an authorized dealer and has to be operated properly to be warranted. This warranty is limited to the original purchaser and does not cover defects on the vehicle. The warranty registration must be completed and returned within 30 days of purchase. This warranty does not cover abnormal wear or damages resulting from an accident, misuse, neglect or if the product has been altered in any way or aspect. Defective parts must be returned by an authorised dealer, who will have received a prior return authorization from Adapt Solutions #1 Ltd. All warranted parts shall be sent to Adapt Solutions #1 Ltd prepaid. Parts will be replaced or repaired at Adapt Solutions #1 Ltd.

Date installed: _____

Dealer: _____

Dealer signature: _____

TO KEEP FOR YOUR RECORDS

(Place serial number here)

WE ARE THE SOLUTION

Please contact your dealer for immediate assistance

Dealer:

TF 1.866.641.0419 | tech@adaptsolutions.com



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | Fax 418.889.9838 | P 418.889.0419

www.adaptsolutions.com

