

Evergreen Surfaces

Evergreen has a wide range of surfaces to meet the needs of clients from high risk to clients with existing pressure injuries (pressure injury stages I – IV). They are designed to prevent and manage skin injuries using alternating pressure therapy, low air loss and 'Auto-Detect' technologies. Our surfaces are designed from a clinician's perspective – making them intuitive to use and suitable for all patient care environments.

Alternating Pressure Therapy

1-in-3 air cells alternately inflate and deflate to mimic the movement of the body, enhance pressure redistribution and promote blood flow to high risk tissues allowing greater oxygenation and healing.

Low Air Loss Therapy

Constant flow of air through the surface to maintain constant temperature and reduce humidity. Promotion of skin integrity through reduced moisture build-up.

'Auto-Detect'

Auto-Detect' technology automatically calibrates optimal pressure based on bodyweight and surface area. This improves the efficiency of set-up and reduces the need for on-going adjustments.

Compare Our Surfaces

| | Palm | OpalAir | OpalAir Turn |
|---|-------------------------------|--------------------------------------|------------------------------|
| Stage | I – III | I – IV | I – IV |
| Low Air Loss | Micro | True | True |
| Size Options | 36" (extend with bolsters) | 36" or 42" (extend with bolsters) | 36" or 42" |
| Weight Cap | 350lbs | 600lbs (36") 800lbs (42") | 600lbs (36") 800lbs (42") |
| 'SmoothCell' for stability and comfort | ✓ | ✓ | ✓ |
| 'StaticPillow' motionless head section | ✓ | ✓ | ✓ |
| Quietest in class low noise pump | ✓ | ✓ | ✓ |
| Power Outage Coverage | ✓ | ✓ | ✓ |
| 'EasyCover' easy-to-remove, machine wash | ✓ | ✓ | ✓ |
| 'Auto-Fowler' cells inflate when be head raised | ✗ | ✓ | ✓ |
| Cool/Warm Air function | ✗ | ✓ | ✓ |
| Lateral Rotation Therapy | ✗ | ✗ | ✓ |
| MSRP | Contact your Care Provider | Contact your Care Provider | Contact your Care Provider |